

**Product Name :**  
Pharma Automatic Tablet and Capsule Dissolution Tester

**Product Code :**  
NLAB-PHARMACYLABM11001



**Description :**

Pharma Automatic Tablet and Capsule Dissolution Tester

**Technical Specification :**

It has four vessels and three poles ,arranged on one line ;

The head part of the instrument can be artificially turned over smoothly and flexibly .

The basket units and paddle units are made of imported stainless steel (ss316L)

The magnetism-pumped circular water current can heat the system evenly ,and the bath liquid can achieve an equal temperature.

A MPU is used in dissolution tester to automatically control the temperature,rotational speed and time ,and the parameters can be preset ant any time .

Temperature range: ambient to +39.9°C

Speed range: (20—200)rpm

Speed accuracy:  $\pm 2$ rpm

Temperature Stability:  $\pm 0.3^{\circ}\text{C}$

Automatic time range: (5-900)min

Time accuracy:  $\pm 0.5$ min

Paddle radial run-out:  $\pm 0.5$ min

Basket radial run-ou:  $\pm 0.1$ min

Power: 220v/50HZ/1000W or 110v/60HZ/1000W

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003