

Product Name :
Pharma Automatic Tablet and Capsule Dissolution Tester

Product Code :
NLAB-PHARMACYLABM11001



Description :

Pharma Automatic Tablet and Capsule Dissolution Tester

Technical Specification :

It has four vessels and three poles ,arranged on one line ;

The head part of the instrument can be artificially turned over smoothly and flexibly .

The basket units and paddle units are made of imported stainless steel (ss316L)

The magnetism-pumped circular water current can heat the system evenly ,and the bath liquid can achieve an equal temperature.

A MPU is used in dissolution tester to automatically control the temperature,rotational speed and time ,and the parameters can be preset ant any time .

Temperature range: ambient to +39.9°C

Speed range: (20—200)rpm

Speed accuracy: ±2rpm

Temperature Stability: ±0.3°C

Automatic time range: (5-900)min

Time accuracy: ± 0.5 min

Paddle radial run-out: ± 0.5 min

Basket radial run-ou: ± 0.1 min

Power: 220v/50HZ/1000W or 110v/60HZ/1000W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003