

Product Name :
Pharma Automatic Tablet Hardness Tester

Product Code :
NLAB-PHARMACYLABM11003



Description :

Pharma Automatic Tablet Hardness Tester

Technical Specification :

The tablet is loaded and pressed artificially and pressed automatically.

There is high accuracy press sensor and hardness data is displayed on LED.

It tests tablet hardness continuously. The biggest number is 100 pieces of a gang.

Time of testing will be saved by directly measuring the diameter .

Test data is displayed and locked automatically .so can test be reset and cycled .

It can deal with and print out such data as : hardness of each tablet ,the maximum value ,the minimal value ,the average value and the number of tablets.

The parameters can be preset at any time.

Hardness range: (2 – 199. 9) N

Hardness accuracy: Max. * (+ 0.05 %)

Moving distance of probe: 3-40mm

Test period: 20 S/P

Max number in each group: 100P

Power: 220V /50HZ/60W or 110V /60Hz/60W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003