

Product Name :
Pharma Dissolution Testing Machine 6 Cups

Product Code :
NLAB-PHARMACYLABM11002



Description :

Pharma Dissolution Testing Machine (6 Cups)

Technical Specification :

Temperature range: ambient to +39.9°C

Speed range: (20—200)rpm

Speed accuracy: ± 2 rpm

Temperature Stability: $\pm 0.3^\circ\text{C}$

Preset timing: six point (1-999)min

Time accuracy: ± 0.5 min

Paddle radial run-out: ± 0.5 min

Basket radial run-ou: ± 0.1 min

Power: 220v/50HZ/1000W or 110v/60HZ/1000W

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003