

**Product Name :**  
Pharma Gas Filling Machine

**Product Code :**  
NLAB-PHARMACYLABM13002



**Description :**

Pharma Gas Filling Machine

**Technical Specification :**

For valves: 1 " (25.4 mm)

Contents of the filling volume (ml): 20-450

Contents of the input connector: 20 mm diameter stainless steel pagoda connector

Propellant filling volume (ml): 20-450

Propellant input connector: M22 \* 1.5 Anti-teeth (for example, with standard 15KG, 50KG liquefied petroleum gas cylinder butt)

Filling accuracy:  $\pm 1\%$

Aerosol cans for height (mm): 70-330

For aerosol can diameter Diameter (mm):  $\varnothing 35$  to  $\varnothing 70$

Working pressure (MPa): 0.5-1.0

Maximum Air Consumption (m<sup>3</sup>/min): 1

Production capacity (cans / hour): 100 ml per tank filling of liquid propellant calculations, 1200 cans / hour

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003