

Product Name :
Pharma Semi-Automatic Aerosol Can Filling Machine

Product Code :
NLAB-PHARMACYLABM13001



Description :

Pharma Semi-Automatic Aerosol Can Filling Machine

Technical Specification :

For valves: 1 " (25.4 mm)

Contents of the filling volume (ml): 20-450

Contents of the input connector: 20 mm diameter stainless steel pagoda connector

Propellant filling volume (ml): 20-450

Propellant input connector: M22 * 1.5 Anti-teeth (for example, with standard 15KG, 50KG liquefied petroleum gas cylinder butt)

Filling accuracy: $\pm 1\%$

Aerosol cans for height (mm): 70-330

For aerosol can diameter Diameter (mm): $\varnothing 35$ to $\varnothing 70$

Working pressure (MPa): 0.5-1.0

Maximum Air Consumption (m³/min): 1

Production capacity (cans / hour): 100 ml per tank filling of liquid propellant calculations, 1200 cans / hour

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003