

Product Name :
Pharma Speed Tablet Blister Packing Machine

Product Code :
NLAB-PHARMACYLABM15001



Description :

Pharma Speed Tablet Blister Packing Machine

Technical Specification :

Cutting Frequency: Al-plastic:25-50times/min

Production Capacity: Al/plastic:120thousand Grains/hour;

Max Forming Depth: Al/plastic:20mm (customized forming depth)

Travel Range: 30-120 mm (made to order)

Total Power: 380V 50HZ 8.5KW

PVC Rigid Sheet: 0.25-0.5x260 mm

Dialysis Paper: 50-100gx260 mm

PTP Aluminum Foil: 0.02-0.035 x260 mm

Cooling: Condensate(self)or cold circulating water(optional)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003