

**Product Name :**  
Pharmaceutical Blister Packing Machine

**Product Code :**  
NLAB-PHARMACYLABM15003



**Description :**

Pharmaceutical Blister Packing Machine

**Technical Specification :**

Capacity: 20,000-50,000grain/hour

Punching frequency: 30times/min/plate each punching

Plate size: 90x80mm(Designed according user's requirement)

Packing material: PVC hard piece:width65-90mm,thickness0.25-0.30mm,bore diameter of rolling shaft60-70mm

PTP Aluminum foil: width65-70mm,thickness0.02-0.025mm,bore diameter of rolling shaft60-70mm

Power consumption: 5kw,220v50-60hz,1-pole 5-line

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

