

**Product Name :**  
Photo Diode AS Measurement Transducer Trainer

**Product Code :**  
NLAB-ELECTRONICSAB110020



**Description :**

Photo Diode AS Measurement Transducer Trainer

**Technical Specification :**

**SCOPE OF LEARNING:**

Study of Linear Displacement Measurement Using Photo Diode

**TECHNICAL SPECIFICATIONS:**

**Power Supplies:**

DC Dual Power Supply IC Regulated +12V DC, 150mA.

DC Power Supply IC Regulated +5V DC, 150mA.

DC Power Supply IC Regulated 0-12V DC, 150mA.

Operated on Mains power 230V, 50Hz +10%

**Digital Meters:**

Voltmeter 200mV (Distance in mV)

Components are mounted on the panels are:

Resistors

Capacitors

IC 741

Variable Resistor For Gain and Zero Setting

Differential Amplifier With Feedback

Lead Compensation System

O/P Provided On Test Points for Monitoring & Controlling.

Photo Diode Sensor Mounted on External Base with Output Carrying out through 2mm Connector.

#### SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003