

Product Name :
Photo Stability Test Chamber

Product Code :
NLAB-LABORATORYEM32005



Description :

Photo Stability Test Chamber

Technical Specification :

Photo Stability Test Chamber These high-quality photo stability test chamber will have a diffuser module with protection coating to offer a UV light energy, easily manageable from 200 watt hour per sq m.

These superior technological arrangements enable a direct evaluation between the drug substance and drug product to give a fair result.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003