

Product Name :
Physics Lab Ammeter Laboratory Type

Product Code :
NLAB-EDUSCHOOLEQ158059



Description :

Physics Lab Ammeter Laboratory Type

Technical Specification :

Range: 0-50mA-500mA-5A

Compact size

On base square shape with terminal

Ammeter Type: DC

Easy operation and stable work, Binding post connection

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003