

Product Name :
Physics Lab Deflection Magnetometer

Product Code :
NLAB-EDUSCHOOLEQ158052



Description :

Physics Lab Deflection Magnetometer

Technical Specification :

It is used for the measurement of the Earth's magnetic field, in military for measuring the submarines and also detects magnetic anomalies of various types.

Compass box fitted on polished wooden base made of non-wrapping wood with meter scale alongside its length.

Compass box is removable.

Magnetometer is used for determination of magnetization of a magnetic material, the strength and also the direction of the magnetic field at a point in space.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003