

**Product Name :**  
Physics Lab Electromagnet

**Product Code :**  
NLAB-EDUSCHOOLEQ158054



**Description :**

Physics Lab Electromagnet

**Technical Specification :**

Electromagnet U shape has two coils of copper and inserted in U shape iron rod.

Comprising iron U-shaped core of circular section with armature and carrying hook.

Wound with copper wire soft iron U-shaped core

Core size : 100 mm long and 10 mm dia. (Min)

Armatures on both the arms wound with enameled copper wire and are joined in series.

Operates on 3-6 Volts D.C.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

