

Product Name :
Physics Lab Fluid Pressure Apparatus

Product Code :
NLAB-EDUSCHOOLEQ158017



Description :

Physics Lab Fluid Pressure Apparatus

Technical Specification :

Features two brass cylinders (19mm and 38mm diameter) fitted with pistons and connected by a brass tube.

The pistons are fitted with circular discs for loading masses (masses not included).

The apparatus is fixed to a sturdy molded plastic base with rubber stoppers underneath to prevent the apparatus from sliding or moving while in use.

To demonstrate the transmissibility of fluid pressure.

Made of two metal cylinders, Cylinders dia.: 15 mm and 30 mm respectively

Cylinders fitted with piston and connected by a tube.

Complete system would be fitted on a wooden base.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments