

**Product Name :**  
Physics Lab Hydrostatic Balance

**Product Code :**  
NLAB-EDUSCHOOLEQ158020



**Description :**

Physics Lab Hydrostatic Balance

**Technical Specification :**

Capacity: 0 - 200 gm

Sensitivity: 10 mg (Min)

Pan diameter: 75 mm (Min)

Fine Marble base with spirit level and leveling screw in glass fitted case, complete with weight set

Finished in durable black enamel with chrome plated fitting and mounted on nicely polished heavy wooden base provided with leveling screws.

A third small pan with hook is supplied to use it as hydrostatic balance for Archimedes's principle.

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## **Naugralabequipments**

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)