

Product Name :
Physics Lab Milli Ammeter AC (Analog)

Product Code :
NLAB-EDUSCHOOLEQ158064



Description :

Physics Lab Milli Ammeter AC (Analog)

Technical Specification :

2 Ranges: 10/100mA

Accuracy: $\pm 2.5\%$ F.S. Deflection Angle: 90°

Scale Length: 65mm

Terminal: 4mm socket captive head, suitable for wire or pin-type test lead.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003