

**Product Name :**  
Physics Lab Milli Ammeter

**Product Code :**  
NLAB-EDUSCHOOLEQ158061



**Description :**

Physics Lab Milli Ammeter

**Technical Specification :**

2 Ranges: 50/500mA, Accuracy:  $\pm 2.0\%$  F.S.

Deflection Angle:  $90^\circ$

Scale Length: 65mm

Terminal: 4mm socket captive head, suitable for wire or pin-type test lead.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003