

**Product Name :**  
Physics Lab Rheostats

**Product Code :**  
NLAB-EDUSCHOOLEQ158067



**Description :**

Physics Lab Rheostats

**Technical Specification :**

Good linearity, Sliding contact of coppers graphite, Small and light

Resistance tolerance:  $\pm 10\%$

Open slide wire type with a variety of resistance and current carrying capacity.

Copper nickel alloys wire is oxidized to provide perfect insulation.

Earthing resistance:  $< 0.1 \text{ Ohm}$ , Rated range: 5 Ohm, 3A

Max. working voltage: 380VAC, 400VDC

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

