

Product Name :
Physics Lab Sonometer

Product Code :
NLAB-EDUSCHOOLEQ158046



Description :

Physics Lab Sonometer

Technical Specification :

Comprising a hollow wooden box with meter rules divided in mm.

Three wires of brass & steel attached to fixed bridge & tensioned by fine adjustment pegs.

Two movable bridges, pulley & fixed screw to enable a third wire to fixed for tensioning by masses.

Would complete with two bridges.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003