

**Product Name :**  
Physics Lab Thermometer

**Product Code :**  
NLAB-EDUSCHOOLEQ158030



**Description :**

Physics Lab Thermometer

**Technical Specification :**

Glass Thermometer

Material : Glass

Used in laboratory

Graduated

Temperature Range : upto 10°C to 150°C, 0°F to 400°F, rooms 0°R to 80°R Min

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003