

Product Name :
Physics Lab Tool Box

Product Code :
NLAB-EDUSCHOOLEQ158075



Description :

Physics Lab Tool Box

Technical Specification :

The multi-purpose tool kit consisting of the following tools :-

Long nose pliers, Diagonal cutting pliers, Fine tip straight tweezers, Soldering aid tools (3pcs), Adjustable wrench, Ceramic soldering iron, Screwdriver handle, Pro-soft s/d 5.0x75mm and #1x75mm, De-soldering wick, Crimping tool (metric), De-soldering pump, IC extractor and IC insert tool, Solder core 63% SN, Prong holder and Parts tube, Heat sink and Bits pack (8pcs), Zipper carrying bag

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003