

Product Name :
Physics Lab Tuning Fork

Product Code :
NLAB-EDUSCHOOLEQ158044



Description :

Physics Lab Tuning Fork

Technical Specification :

Set of 8

Plain shank

Material : Steel

Packed in wooden box

Frequencies : 256Hz, 288Hz, 320Hz, 341.3 Hz, 384Hz, 426.6Hz, 480Hz & 512Hz

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003