

Product Name :
Physics Lab Vacuum Pump

Product Code :
NLAB-EDUSCHOOLEQ158025



Description :

Physics Lab Vacuum Pump

Technical Specification :

V- Belt driven and mounted on base plate with belt guard.

Vacuum pump with motor 220volt AC single phase.

Vacuum displacement 32 litre / Min.

Compact Design

Low Noise Level

Maintenance Free

Easy to handle

Continuous Duty

Totally Dry Compressed Air

Low Running Cost

Pump Speed (m3/h): 3.5 (Min)

Ultimate Pressure (MPa): 0.09 (Min)

Rotate Speed: 1300 rpm (Min)

Working Voltage: 220 V / 50 Hz (Min)

Motor Power: 0.25 KW (Min)

Inlet Caliber: Dia. 7 mm (Min)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003