


Product Name :
Physics Lab Vernier Caliper SS

Product Code :
NLAB-EDUSCHOOLEQ158006



Description :

Physics Lab Vernier Caliper SS

Technical Specification :

Plated steel, Columbus type caliper enabling the user to read either inside or depth measurement.

Material: Stainless Steel

Measuring range: 0-300mm

Reading: 0.1mm or 0.05mm

Can be measured OD, ID, depth and steps, Would have a lock screw for holding the sliding jaw position.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003