

Product Name :
Physics Lab Voltmeter AC (Analog)

Product Code :
NLAB-EDUSCHOOLEQ158065



Description :

Physics Lab Voltmeter AC (Analog)

Technical Specification :

Range: 0-15V-150V

Compact size

A voltmeter, also known as a voltage meter, is an instrument used for measuring the potential difference, or voltage, between two points in an electrical or electronic circuit.

It is generally used in all types of Laboratories such as School, college, university, Research etc.

Easy operation and stable work, Binding post connection

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

