

Product Name :
Physics Lab Voltmeter

Product Code :
NLAB-EDUSCHOOLEQ158062



Description :

Physics Lab Voltmeter

Technical Specification :

Range: 0-3V-30V-300V

Compact size

It is generally used in all types of Laboratories such as School, college, university, Research etc.

DC mounted in beautiful Bakelite / acrylic stand with accuracy +0% F.S.D.

Voltmeter Type: DC

Easy operation and stable work

Binding post connection

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments