


**Product Name :**  
PIC16F Training Lab

**Product Code :**  
NLAB-ENGINEERINGLB49004



**Description :**

PIC16F Training Lab

**Technical Specification :**

PIC16F Training Lab which uses Microchip's PIC16F887 is a 8-bit microcontroller to implement various I/O control experiments. This controller contains most of the powerful functions in modern MCUs nowadays. Moreover, it can be used for automation, motor control, device measurement, and mechanical controls...etc. PIC16F is popular and well-known due to its economic cost, wide applicability, high accessibility and reliable stability.

contains several peripheral devices from basic LED to advanced capacitive sensing module. The integration of these devices enable end users to create and implement different kind of control experiments. Together with our user friendly experiment manual, end users can learn the control of PIC MCU more conveniently and efficiently

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

