


**Product Name :**  
Piezo Electric Transducer Trainer

**Product Code :**  
NLAB-ELECTRONICSAB110017



**Description :**

Piezo Electric Transducer Trainer

**Technical Specification :**

**SCOPE OF LEARNING:**

Study of Piezo Electric Transducer

**TECHNICAL SPECIFICATIONS:**

**Digital Meters:**

Voltmeter 200mV DC.

**Power Supplies:**

DC Supply IC Regulated +12V DC, 150mA.

DC Supply IC Regulated +5V DC, 150mA.

Operated on Mains power 230V, 50Hz +10%

Components are mounted on the panels are:

Variable Resistor (Presets)

741 IC

Piezo Electric Sensor.

**SALIENT FEATURES:**

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Instruction manual.

Protection Covers.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## **Naugralabequipments**

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003