


Product Name :
Pilling Tester Machines

Product Code :
NLAB-TESTINGLABM64001



Description :

Pilling Tester Machines

Technical Specification :

TECHNICAL DATA

Size of wooden boxes (Inside) 225 x 225 x 225 mm

Thickness of cork lining 3 mm

Speed of rotation of boxes 60 rpm

Size of test specimen 125 x 125 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003