

Product Name :
Pipe Bender

Product Code :
NLAB-VOCATIONALTVET16003



Description :

Pipe Bender

Technical Specification :

Steel construction.

6 cars iron shoes for 0.5 - 2 in. diam. pipe.

Makes 90-degree bends.

Top rollers adjust for degree of bend.

Overall dimensions: 22.62L x 21.26W x 6.3H in.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003