

**Product Name :**  
Pipe Bending Machine

**Product Code :**  
NLAB-VOCATIONALTVET16006



**Description :**

Pipe Bending Machine

**Technical Specification :**

Capacity	
250 mm Stroke	
215-680 mm (mm.)	

Bending Frames (2, 2)½, 3, 4
275 Thick 6.5mm Pipe
700mm (203)
500kg (Approx.)

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003