

Product Name :
Pipe Bending Machine

Product Code :
NLAB-VOCATIONALTVET16006



Description :

Pipe Bending Machine

Technical Specification :

Capacity	
250 mm Stroke	
215-685 mm (mm.)	

2. Bending Form (2, 2) 1/2, 3, 4
2. 25 mm thick 6.5 mm Pipe
2. 25 mm (20)
2. 25 mm (Approx.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003