

Product Name :
Pipe Loss Test Bench

Product Code :
NLAB-EDUCATIONTVET14005



Description :

Pipe Loss Test Bench

Technical Specification :

This equipment makes it possible to measure the pressure drop across a pipe

The bench is composed of:

Several test tubes, length 1 m, internal diameters 4.5; 7.7; 10.9; 17.2mm

An artificially roughened 15.2 mm ID tube

Elbows 90 °, 45 °, Y 45 °, T 90 °

Diameter enlargements and reductions

Different valves

Transparent pipe section with Venturi, orifice flowmeter and mounted Pitot probe

Accessories

Digital pressure gauge

1m pressurized water pressure gauge

Total position 5

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003