

**Product Name :**  
Pipe Threading Tool

**Product Code :**  
NLAB-VOCATIONALTVET16001



**Description :**

Pipe Threading Tool

**Technical Specification :**

Pipe Threading Tool

Four Cast iron die heads with alloy steel dies.

Two-piece tubular steel handle.

Compact ratchet head.

The compact ratchet head is suitable for working close to walls.

Packed in blow mould carrying case and display sleeve.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003