

Product Name :
Pipe Welding Clamp

Product Code :
NLAB-VOCATIONALTVET13038



Description :

Pipe Welding Clamp

Technical Specification :

Aluminum alloy frame for high strength and lightweight.

Four self-centering clamps.

No adapters needed.

Compatible with tapping tees.

Fixed positions for 0, 45, 90 degrees

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003