

Product Name :
Pipette Pump

Product Code :
NLAB-EDUSCHOOLEQ126001



Description :

Pipette Pump

Technical Specification :

Polylab presents a fast release pipette pump device for precise pipetting & zip quick emptying.

The knurled thumbwheel can be rotated to draw up or dispense liquids while the releaser may be pushed in for zip-quick emptying.

The moulded chuck has silicon rubber collet inside to hold various pipettes of standard sizes.

It resists acids & alkalies and can be easily disassembled for thorough cleaning.

Packing

22 Pcs

10 Pcs

22 Pcs

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003