

**Product Name :**  
Piston Pump Cutaway

**Product Code :**  
NLAB-AUTOMOTIVETVET21003



**Description :**

Piston Pump Cutaway

**Technical Specification :**

This cutaway model is carefully sectioned for training purposes, professionally painted with different colours to better differentiate the various parts and cross-sections. Many parts have been chromium-plated and galvanized for a longer life.

Automotive Technologies, Cutaways Models, Acrylic Static Cutaways Equipments for Automotive Training Lab

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003