

Product Name : Pitot Tubes	Product Code : NLAB-TECHNICALAB26010
Description :	
Pitot Tubes	
Technical Specification :	
<p>150mm Pitot tube 150mm traverse Pitot tube supplied with connectors and 10m tubing. 300mm Pitot tube 300mm traverse Pitot tube supplied with connectors and 10m tubing. 450mm Pitot tube 450mm traverse Pitot tube supplied with connectors and 10m tubing.</p> <p>Tubes are in stainless steel and mounted on a supporting body with scale. They are supplied with a watertight gland for installation below water level. In order to measure velocity, the Pitot tubes must be connected to a manometer.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003