

Product Name :
Plancks Constant Determination Box Demonstration Kit

Product Code :
NLAB-SCHOOLABS17009



Description :

Plancks Constant Determination Box Demonstration Kit

Technical Specification :

Determination of Planck's constant Using the working principle of LED lamps to determine the wavelength
Covers a wavelength of 470 – 940 nm (7 ranges).

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003