

**Product Name :**  
Planetary Mixer Table Top

**Product Code :**  
NLAB-VOCATIONALTVET27009



**Description :**

Planetary Mixer Table Top

**Technical Specification :**

Body Material: Stainless Steel

Automation Grade: Semi-Automatic

Capacity: 7 Liter

Motor Power: 1 Hp

Power Source: Electric

Motor Speed: 2000 RPM

Machine Type: Table Top Mixer

Voltage: 220 Volt

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003