

Product Name : Planetary Mixer	Product Code : NLAB-VOCATIONALTVET27002
Description :	
Planetary Mixer	
Technical Specification :	
<p>The mixer has two paddles and one scraper that help process the food quickly and minces efficiently.</p> <p>The mixer is equipped with three types of a tools namely flat beater, whisk, and J hook.</p> <p>The role of the flat beater includes mix base of the products, whisk allow mixing of the air in the cream base and J Hook mixes the hard dough.</p> <p>Their three-speed leveler provided that is easy to shift.</p> <p>Features:</p>	

Suitable for producing liquids, ointments, hard food, powder and other material

Can operate using steam, electricity, and ga

Hydraulic lifting, stable and reliable

Specifications:-

Body Material: Stainless Steel

Automation Grade: Automatic

Color: White

Phase: Single Phase

Motor Power: 0.5 Hp

Max. Loading Capacity: 3 Litres

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003