

**Product Name :**  
Plasma Welding Machine

**Product Code :**  
MU144



**Description :**

Plasma Welding Machine

**Technical Specification :**

Plasma Welding Machine  
With Coolant recirculatory,  
Cooling capacity 1.75gal.  
Auto/Manual Operation,  
Pilot Arc.  
220V/60HZ Single Phase.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003