

Product Name :
Plastic Box

Product Code :
NLAB-EDUSCHOOLEQ46021



Description :

Plastic Box

Technical Specification :

Plastic Box is used as a storage container.

It can be either used to store samples of different materials as well as components like wires, screws, tips or small tools etc.

It is made of clear plastic with arrangement of bottom locking the lid so that the materials in the box remain secure inside it.

These transparent storage boxes are available in different sizes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003