

**Product Name :**  
Plastic Limit Set

**Product Code :**  
NLAB-VOCATIONALTVET18019



**Description :**

Plastic Limit Set

**Technical Specification :**

Glass plate 20 cm x 15 cm having ground ends and one side frosted.

Brass or stainless steel rod 3 mm Dia x 150 mm long.

Flexible Spatula.

Set of 6 moisture containers.

Porcelain basin 150 mm Dia.

Plastic wash bottle 500ml.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

