

Product Name :
Plastic Limit

Product Code :
NLAB-TESTINGLABM109001



Description :

Plastic Limit

Technical Specification :

Shrinkage limit is the maximum water content at which a reduction in water content does not cause an appreciable reduction in volume of the Soil Mass. At shrinkage limit, on further reduction in water, air enters in to the voids of soils and thus keeps the volume constant. The apparatus can be used to determine shrinkage limit and to calculate other shrinkage ratio, shrinkage index and volumetric shrinkage.

Specification :

Set consists of one each :-

- i) Porcelain evaporating dish,
- ii) Shrinkage dish,
- iii) Glass cup,

iv) Perspex plate with three metal prongs,

v) Flexible spatula 100mm,

vi) Glass cylinder,

vii) 25ml x 1ml, supplied without mercury

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003