

Product Name :
Plasticity Penetrometer for Liquidity Limit

Product Code :
NLAB-EDUCATIONTVET15055



Description :

Plasticity Penetrometer for Liquidity Limit

Technical Specification :

Solid base with column on which is fixed a cone penetrometer with a smooth surface at 30 °, and comprising an indicator dial of d = 150mm, graduated in divisions of 0.1mm.

The height of the cone is easily adjusted according to the size of the sample. Reading the penetration is simple and replacing the cone is quick and easy.

Bucket for penetration test of d = 55mm and h = 40mm. Edge parallel to the base.

Gauge to test the tip condition of the cone

Mirror for easy adjustment

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

