

**Product Name :**  
Plate Bearing Test

**Product Code :**  
NLAB-VOCATIONALTVET18012



**Description :**

Plate Bearing Test

**Technical Specification :**

Hydraulic jack with fixed pumping unit capacity: 0-500 kN. x 0.5 kN.

Pressure gauge and flexible metal pipe 5m long - 1No.

1 Special ball and socket arrangement between the jack and the bearing plate.

16 Extension rod 12mm. dia x 25cm long.

4 Magnetic base with a female thread on top for holding extension rod.

4 Top end plate, 50mm, dia with a male thread for fitting onto the extension rods and positioning the dial gauge plunger.

2 Column 15cm dia x 25cm long with flanges complete with four bolts and nuts.

1 Column 15cm x 50cm long with flanges complete with four bolts and nuts.

1 Bridge support for welded steel angle construction, 5m, span and stands approximately 30cm, high, fitted with two quick release clamps, dial bracket.

Plane M.S (Mild Steel) Plate 60 cms x 60cms sq. x 25mm thick.

Plane M.S.Plate 45cms x 45cm.sq. x 25mm thick.

Plane M.S.Plate 30cms x 30cm.sq. x 25mm thick.

4 Dial Gauge 0.01mm x 25mm.

Accessories:-

Plane Mild Steel (M.S) Plate 75cm x 25mm thick.

Plane M.S. Plate 50cm x 25mm thick.

Grooved M.S. Plate:-

60cm x 60cm x 25mm thick,

45cm x 45cm x 25mm thick,

30cm x 30cm x 25mm thick,

75cm x 75cm x 25mm thick.

Flexible rubber pipe 1m. Length for the hydraulic jack in the place of metal pipe

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003