


Product Name :
Plate Load Bearing Test Apparatus

Product Code :
NLAB-SOILELAB140006



Description :

Plate Load Bearing Test Apparatus

Technical Specification :

Plate Loading Test Set consists of a 200 kN capacity piston assembly and hydraulic hand pump, pressure transducer, LPI battery operated digital readout unit, 2.4 m long datum bar, 3 pcs. 25 mm travel x 0.01 mm digital dial gauges with dial supports, 300 mm and 450 mm dia. loading plates.

Plate Loading Test Set Set consists of a 500 kN capacity piston assembly and hydraulic hand pump, pressure transducer, LPI battery operated digital readout unit, 2.4 m long datum bar, 3 pcs. 25 mm travel x 0.01 mm digital dial gauges with dial supports, 600 mm and 760 mm dia. loading plates.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003