

Product Name :
Plate Rolling Machine

Product Code :
NLAB-MECHANICALAB260010



Description :

Plate Rolling Machine

Technical Specification :

3 Roll PYRAMID Type outline with twofold pre-squeezing game plan

Top roll is repaired and base moves and down.

Base moves up-down through water powered chambers.

Two base rolls are force driven via planetary apparatus box with electric powered engine.

Twisting rolls are made of combination steel.

Dropper end assurance by pressure driven barrels.

Over-burden prevention

Assembled – in cone bowing connection.

Global quality Hydraulic and Electrical.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003