

Product Name : Platelet Agitator	Product Code : NLABMEDICAL-M3200003
Description :	
Platelet Agitator India, Platelet Agitator manufacturer India, suppliers, Medical Platelet Agitator India, Hospital Platelet Agitator India	
Technical Specification :	
Platelet agitator specially designed to maintain, platelets in continues motion at a fixed speed of 70 RPM. Maximum 36 standard Platelet bags can be accommodate at a time.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003